

# Signal Acquisition Team

## Spring 2025 Semester Project

### Real-Time Data Validation Pipeline

Longhorn Neurotech - Software Division

January 2025

## 1 Project Overview

Your mission is to implement a **real-time data validation layer** that catches hardware and signal quality issues *before* data enters the processing pipeline. This is the first line of defense against poor data quality.

### 1.1 Deliverables

- `data_validation.py` - New module with all validation functions
- Integration into `eeg_processor.py` at specified lines
- No modification of existing filtering code
- Documentation of validation metrics

## 2 Technical Background

### 2.1 Why Data Validation Matters

In real-time BCI systems, bad data causes:

- False classifications (user thinks left, system reads right)
- Model training on garbage data
- Inconsistent performance across sessions

Your validation layer detects:

1. **Hardware failures** - disconnected electrodes, broken channels
2. **Signal anomalies** - flatlines, extreme noise, NaN values
3. **Data integrity** - duplicate channels, missing samples

## 2.2 Validation Happens in Two Places

- **Per-sample validation** - Fast checks on every new data chunk (Lines 151-152 in eeg\_processor.py)
- **Buffer validation** - Comprehensive checks on accumulated data (Line 177 in eeg\_processor.py)

## 3 Component 1: NaN and Inf Detection

### 3.1 What It Is

NaN (Not a Number) and Inf (Infinity) values appear when:

- Division by zero in calculations
- Hardware read errors
- Buffer overflow

### 3.2 Mathematical Definition

For data matrix  $X \in R^{C \times T}$  where  $C$  = channels,  $T$  = time points:

$$\text{Valid}(X) = \begin{cases} \text{True} & \text{if } \forall i, j : X_{ij} \in R \text{ and } |X_{ij}| < \infty \\ \text{False} & \text{otherwise} \end{cases} \quad (1)$$

### 3.3 Implementation

Listing 1: Function: detect\_nan\_inf

```
1 def detect_nan_inf(data):
2     """
3     Detect NaN or Inf values in EEG data.
4
5     INPUT:
6         data: np.ndarray, shape (channels, samples)
7             Raw EEG data from board
8
9     OUTPUT:
10        is_valid: bool
11                True if no NaN/Inf found
12        nan_channels: list
13                    Indices of channels containing NaN/Inf
14
15    EXAMPLE:
16        >>> data = np.array([[1.0, 2.0, np.nan], [3.0, 4.0, 5.0]])
17        >>> is_valid, bad_ch = detect_nan_inf(data)
18        >>> print(is_valid) # False
19        >>> print(bad_ch)   # [0]
20    """
21    # YOUR IMPLEMENTATION HERE
22    # Hint: Use np.isnan() and np.isinf()
23    pass
```

**Your Task:**

1. Check each channel for NaN values using `np.isnan()`
2. Check each channel for Inf values using `np.isinf()`
3. Return list of bad channel indices
4. Set `is_valid = False` if any bad values found

## 4 Component 2: Flatline Detection

### 4.1 What It Is

A **flatline** occurs when a channel outputs constant values, indicating:

- Disconnected electrode
- Hardware malfunction
- Saturated amplifier

### 4.2 Detection Method

A channel is flatline if its standard deviation is below threshold:

$$\text{Flatline}(x) = \begin{cases} \text{True} & \text{if } \sigma(x) < \theta_{\text{flat}} \\ \text{False} & \text{otherwise} \end{cases} \quad (2)$$

where  $\sigma(x)$  is standard deviation and  $\theta_{\text{flat}} = 0.1 \mu\text{V}$  (typical threshold).

### 4.3 Implementation

Listing 2: Function: `detect_flatline`

```

1 def detect_flatline(data, threshold=0.1):
2     """
3     Detect flatline channels (constant/near-constant signal).
4
5     INPUT:
6         data: np.ndarray, shape (channels, samples)
7             EEG data segment
8         threshold: float
9             Standard deviation threshold in microvolts
10            Default: 0.1 uV
11
12     OUTPUT:
13         flatline_channels: list
14             Indices of flatline channels
15         channel_stds: np.ndarray
16             Standard deviation of each channel
17
18     REASONING:
19         Healthy EEG has std > 1 uV typically
20         Flatline means electrode disconnected or broken

```

```

21     """
22     # YOUR IMPLEMENTATION HERE
23     # Hint: Calculate std per channel, compare to threshold
24     pass

```

## 5 Component 3: Extreme Noise Detection

### 5.1 What It Is

**Extreme noise** appears as:

- Muscle artifacts (EMG contamination)
- Movement artifacts
- Electrical interference spikes

### 5.2 Detection Method

Use  $z$ -score to detect outliers:

$$z_{ij} = \frac{X_{ij} - \mu_i}{\sigma_i} \quad (3)$$

Spike detected if  $|z_{ij}| > \theta_z$  where  $\theta_z = 5$  (typical threshold).

### 5.3 Implementation

Listing 3: Function: detect\_extreme\_noise

```

1  def detect_extreme_noise(data, z_threshold=5.0):
2      """
3      Detect extreme noise/spike artifacts.
4
5      INPUT:
6          data: np.ndarray, shape (channels, samples)
7          z_threshold: float
8                      Z-score threshold for spike detection
9                      Default: 5.0 (5 standard deviations)
10
11      OUTPUT:
12          has_spikes: bool
13                     True if spikes detected in any channel
14          spike_channels: list
15                       Channels containing spikes
16          spike_percentage: float
17                          Percentage of samples that are spikes
18
19      DETECTION LOGIC:
20          For each channel:
21              1. Calculate mean and std
22              2. Compute z-scores: z = (x - mean) / std
23              3. Flag samples where |z| > threshold
24              4. If > 1% of samples are spikes, mark channel bad

```

```

25     """
26     # YOUR IMPLEMENTATION HERE
27     pass

```

## 6 Component 4: Channel Duplication Check

### 6.1 What It Is

Sometimes hardware errors cause:

- Same channel duplicated in multiple positions
- Copy of previous buffer instead of new data

### 6.2 Detection Method

Use correlation between channels:

$$\rho(x_i, x_j) = \frac{\text{cov}(x_i, x_j)}{\sigma_{x_i} \sigma_{x_j}} \quad (4)$$

If  $\rho > 0.99$ , channels are duplicates.

### 6.3 Implementation

Listing 4: Function: detect\_channel\_duplication

```

1 def detect_channel_duplication(data, correlation_threshold=0.99):
2     """
3     Detect if channels are duplicates of each other.
4
5     INPUT:
6         data: np.ndarray, shape (channels, samples)
7         correlation_threshold: float
8             Pearson correlation threshold
9             Default: 0.99
10
11     OUTPUT:
12         duplicate_pairs: list of tuples
13             [(ch_i, ch_j), ...] duplicate channel pairs
14         correlation_matrix: np.ndarray, shape (channels, channels)
15             Full correlation matrix
16
17     USAGE:
18         If duplicates found, data is invalid - hardware issue
19     """
20     # YOUR IMPLEMENTATION HERE
21     # Hint: Use np.corrcoef(data) to get correlation matrix
22     pass

```

## 7 Component 5: Validation Package Class

### 7.1 What It Is

A **ValidationPackage** bundles data with its quality metrics.

Listing 5: Class: ValidationPackage

```
1 class ValidationPackage:
2     """
3     Container for validated EEG data with quality metrics.
4
5     ATTRIBUTES:
6         data: np.ndarray
7             The validated EEG data
8         is_valid: bool
9             Overall validity flag
10        timestamp: float
11            Timestamp when data was collected
12        quality_metrics: dict
13            Dictionary of quality measurements
14            {
15                'has_nan': bool,
16                'has_inf': bool,
17                'flatline_channels': list,
18                'noisy_channels': list,
19                'duplicate_pairs': list,
20                'sampling_rate': float
21            }
22
23        USAGE:
24            package = ValidationPackage(data, timestamp)
25            if package.is_valid:
26                # Safe to process
27                process(package.data)
28            else:
29                # Log the issue
30                log_error(package.quality_metrics)
31        """
32        def __init__(self, data, timestamp):
33            self.data = data
34            self.timestamp = timestamp
35            self.is_valid = True
36            self.quality_metrics = {}
```

## 8 Integration Points

### 8.1 File: eeg\_processor.py

Line 1-10: Add import statement:

```
1 # Add after existing imports
2 from data_validation import (
3     detect_nan_inf,
```

```

4     detect_flatline,
5     detect_extreme_noise,
6     detect_channel_duplication,
7     ValidationPackage
8 )

```

**Line 151-152:** Add per-sample validation in `get_recent_data()`:

```

1 # Current line 151:
2 eeg_data = data[self.eeg_channels, :]
3
4 # INSERT AFTER LINE 151:
5 # Validate new data chunk
6 is_valid, nan_channels = detect_nan_inf(eeg_data)
7 if not is_valid:
8     print(f"Warning: NaN/Inf detected in channels {nan_channels}")
9     # Skip this chunk or use last known good data
10    return self.processed_data_buffer[:, -int(duration * self.
        sampling_rate):]

```

**Line 177:** Add buffer validation before returning:

```

1 # Current line 177 returns model_input
2 # INSERT BEFORE RETURN:
3
4 # Comprehensive validation on buffer
5 flatline_ch = detect_flatline(recent_data)
6 if len(flatline_ch) > 0:
7     print(f"Warning: Flatline detected in channels {flatline_ch}")
8
9 has_spikes, spike_ch, spike_pct = detect_extreme_noise(recent_data)
10 if has_spikes and spike_pct > 5.0:
11     print(f"Warning: Excessive noise ({spike_pct:.1f}%) in channels {
        spike_ch}")
12
13 # Then continue with existing return statement

```

## 9 Testing Your Implementation

### 9.1 Test Cases

Listing 6: Test Script - `test_validation.py`

```

1 import numpy as np
2 from data_validation import *
3
4 def test_nan_detection():
5     """Test NaN detection"""
6     # Create test data with NaN
7     data = np.random.randn(8, 100)
8     data[2, 50] = np.nan
9
10    is_valid, bad_ch = detect_nan_inf(data)
11    assert not is_valid, "Should detect NaN"

```

```

12     assert 2 in bad_ch, "Should identify channel 2"
13     print("    NaN detection works")
14
15 def test_flatline():
16     """Test flatline detection"""
17     data = np.random.randn(8, 100)
18     data[5, :] = 0.0 # Flatline channel 5
19
20     flat_ch, stds = detect_flatline(data)
21     assert 5 in flat_ch, "Should detect flatline in channel 5"
22     print("    Flatline detection works")
23
24 def test_noise():
25     """Test noise detection"""
26     data = np.random.randn(8, 100)
27     data[3, 10:20] = 50.0 # Large spike
28
29     has_spike, spike_ch, pct = detect_extreme_noise(data)
30     assert has_spike, "Should detect spike"
31     print("    Noise detection works")
32
33 if __name__ == "__main__":
34     test_nan_detection()
35     test_flatline()
36     test_noise()
37     print("\nAll tests passed!")

```

## 10 Expected Outcomes

### 10.1 Performance Metrics

- **Validation time:** < 5ms per chunk
- **False positive rate:** < 1%
- **True positive rate:** > 95% on simulated bad data

### 10.2 Documentation Required

1. Comments explaining each validation function
2. Test results showing validation catches errors
3. Performance benchmarks (timing measurements)
4. Integration checklist confirming all insertion points completed



## 11 Grading Rubric

Component	Points	Criteria
NaN/Inf Detection	15	Correctly identifies bad values, returns channel indices
Flatline Detection	15	Uses std threshold, identifies constant channels
Noise Detection	20	Implements z-score method, reports percentage
Duplication Check	15	Calculates correlations, identifies duplicates
Integration	20	Code properly inserted at specified lines
Testing	10	Test script runs and passes all tests
Documentation	5	Clear comments and docstrings
<b>Total</b>	<b>100</b>	

## 12 Timeline

- **Week 1-2:** Implement validation functions
- **Week 3:** Integration into eeg-processor.py
- **Week 4:** Testing with synthetic bad data
- **Week 5:** Hardware testing with real OpenBCI board
- **Week 6:** Documentation and final report

## 13 Resources

- NumPy documentation: <https://numpy.org/doc/>
- EEG artifact detection paper: Delorme et al. (2007) - EEGLAB artifact rejection
- BrainFlow API: <https://brainflow.readthedocs.io/>